

## STRONGBACK 24HD Heavy Duty



## STRONGBACK 24HD Heavy Duty



### UNIQUE

Patented, adjustable ergonomic backrest for unmatched postural support.



### SUPERIOR QUALITY

Impressive weight capacity of 170 kg with a weight of only 16 kg.



### AWARD-WINNING

The lightweight design enables superior manoeuvrability.

Measurements and Weights	
Name	STRONGBACK 24HD Heavy Duty
Rear wheel size	24 in
Seat width	56 cm
Seat depth	48 cm
Seat height	50 cm
Seat angle	12°
Total width	80 cm
Total length	109 cm
Total height	88 cm
Backrest height	45 cm
Backrest angle	curved from -10° bis +24°
Net weight <sup>1</sup>	12.6 kg
Minimum usable weight <sup>2</sup>	14.4 kg
Total weight <sup>3</sup>	16.3 kg
Maximum load	170 kg

<sup>1</sup>Net weight excluding wheels (QR wheels), legrests and anti-tippers

<sup>2</sup>Minimum usable weight including wheels, excluding legrests and anti-tipper

<sup>3</sup>Total weight inclusive wheels, excluding legrests and anti-tippers

