

SEATA

Ergonomic Rollator



2023 Winner!
The Chicago Athenaeum
Museum of Architecture and Design



Sit in Comfort : Walk with Confidence



STRONGBACK
mobility 

Features



UNIQUE

Full-size seat and ergonomic backrest for unmatched comfort.



INNOVATIVE

Our design enables a stable and upright walking position.



AWARD-WINNING

Lightweight and foldable design for easy handling.

Measurements and Weights	
Model	SEATA Rollator
Seat width	48 cm
Seat depth	38 cm
Seat height	50 cm
Seat angle	12°
Front wheel size	23 cm
Rear wheel size	29 cm
Total width	68 cm
Total length	77 cm
Total height	79-99 cm
Backrest height	48 cm
Backrest angle	25°
Size folded	77 x 25 x 80 cm
Total weight	9.3 kg
Maximum load	136 kg

Sit in Comfort : Walk with Confidence



Strongback Mobility Europe GmbH

Erkrather Str. 401, 40231 Düsseldorf • +49 (0) 211 976 33 839

info@strongbackmobility.de • www.strongbackmobility.de